

## PF series MRI

This unique exam light is designed for use in any MRI facility. Made of durable, non-magnetic materials, this mobile light is ideal for imaging centers, radiology departments, and CAT/PET Scan labs and clinics.

### PF MRI FEATURES

- >>> Tested for MRI compatibility
- >>> Powerful, 50W Halogen light source
- >>> Adjustable jointed arm for easy positioning
- >>> Mobile caster base with easy-glide wheels
- >>> Weighted base for added stability



### TECHNICAL DATA

**Power supply:** 230-240V 50 Hz AC or 120V 60 Hz AC  
**Light intensity:** Up to 43,000 Lux/4,000 Footcandles  
**Light source:** Halogen bulb – Dichroic 12V/50W/10°  
**Light color temp:** Up to 3,000 Kelvin  
**Color:** White – RAL 9002

<b>Configurations:</b>	Mobile
<b>Weight:</b>	18.0 kg

Mobile

Art. No.: **PF28016-M**