

Celestial Star™ MRI

The Celestial Star™ MRI brings the performance of a surgical light into the MRI suite. Excellent flexibility and reach combined with high performance and space saving design makes the Celestial Star™ MRI ideal for these facilities. Unlimited positioning and market leading vertical and horizontal reach allow for optimal functionality around the bore. “No tool” service minimizes risk and the compact design limits obstruction of the work area around the MRI equipment.

CELESTIAL STAR™ MRI FEATURES

- >>> Available in ceiling and mobile configurations
- >>> Tested for MRI compatibility
- >>> 730 mm extension arm (optional 990 mm)
- >>> 720° rotation and head-to-toe coverage
- >>> Drift free balance arm with unlimited positioning
- >>> Easily replaceable, long life 3x35W Halogen bulbs



TECHNICAL DATA

Power supply: 230-240V 50 Hz AC or 120V 60 Hz AC
Light intensity: Up to 64,500 Lux/6,000 Footcandles
Light source: Halogen bulb – Dichroic 12V/3x35W/10°
Light color temperature: Up to 4,100 Kelvin
Medical standards: FDA Class II
Color: White – RAL 9002

Configurations:	Ceiling	Mobile
Weight:	23.6 kg	24.3 kg

Ceiling	Art. No.: CS 2050C-MRI
Mobile	Art. No.: CS2050M-MRI
Color Correction Filter	Art. No.: 18676

Note! 3 filters are required. Increases the color temperature to 4,100 Kelvin.

