

## SunExam™

The SunExam™ Series offers a very flexible arm system with installation options for very tall ceilings if needed. The light features a compact lamp head, excellent color temperature, and high light intensity.

### SUNEXAM™ FEATURES

- >>> Large flood pattern, approximately 150 mm diameter
- >>> Easily replaceable, long life 12V/75W Halogen bulb
- >>> 710 mm articulating arm with integral compression spring
- >>> 585 mm extension arm
- >>> IR absorbing, color-correcting glass filter
- >>> Positioning handle mounted on perimeter of lamp head



### TECHNICAL DATA

**Power supply:** 230-240V 50 Hz AC or 120V 60 Hz AC  
**Light intensity:** Up to 42,000 Lux/3,900 Footcandles  
**Light source:** Halogen Bulb 12V/75W  
**Light color temperature:** Up to 4,300 Kelvin  
**Color:** White – RAL 9002

Configurations:	Ceiling	Wall	Mobile
<b>Weight:</b>	21.6 kg	8.7 kg	33.4 kg

Ceiling	Art. No.: <b>28502</b>
<b>Wall</b>	Art. No.: <b>28500</b>
Mobile	Art. No.: <b>28501</b>

