

SunSpot™

Designed for applications where high intensity and very focused lighting is required, the SunSpot™ offers a full range of options including extended ceiling mounts, dimmer functionality, and adjustable spot size. A fiber optic light source is available as an option as well for attachment of fiber optic headlights or instruments.



SUNSPOT™ FEATURES

- >>> Continuously adjustable spot from 75 to 405 mm diameter
- >>> Sterilizable handle
- >>> Mechanical iris dimmer provides constant purity of light without affecting color temperature
- >>> Fiber optic option provides a four position turret for connection to fiber optic head lamps or instruments (ACMI, Storz, Olympus, Wolf)
- >>> IR absorbing, color-correcting glass filter



TECHNICAL DATA

Power supply: 230-240V 50 Hz AC or 120V 60 Hz AC
Light intensity: Up to 21,500 Lux/2,000 Footcandles
Light source: Halogen bulb 20V/150W
Light color temperature: Up to 3,200 Kelvin
Color: White – RAL 9002

Configurations:	Ceiling	Wall	Mobile	Combo	FO
Weight:	26.7 kg	10.7 kg	31.7 kg	42.2 kg	0.5 kg

Ceiling	Art. No.: 28510
Wall	Art. No.: 28508
Mobile	Art. No.: 28509
Ceiling with Fiber Optic	Art. No.: 28514
Wall with Fiber Optic	Art. No.: 28512
Mobile with Fiber Optic	Art. No.: 28513
Combo Ceiling SunSurgery/SunSpot	Art. No.: 28516
Combo Ceiling SunSurgery/SunSpot with Fiber Optic	Art. No.: 28517

