

## SunSurgery™

SunSurgery™ is ideal for surgical applications where shadow management is crucial. This product offers excellent shadow control, superior depth of field and a higher color temperature for tissue clarity.

The modular arm system offers extension arms for tall ceilings, a track system and the option of dual arms in combination with the SunSpot™ product.



## SUNSURGERY™ FEATURES

- >>> Available in ceiling, dual ceiling, combo ceiling, wall, mobile and battery back up configurations
- >>> Prefocused 150 mm at 1.0 meter
- >>> IR absorbing color-correcting glass cylinder
- >>> Easily replaceable, long life 90W Halogen bulb
- >>> 710 mm articulating arm +/-40° vertical adjustment
- >>> 585 mm extension arm
- >>> Optional ceiling track offers 1500 mm of extended reach
- >>> Sterilizable handle



Mobile

### TECHNICAL DATA

**Power supply:** 230-240V 50 Hz AC or 120V 60 Hz AC

**Light intensity:** Up to 72.000 Lux/6,600 Footcandles

**Light source:** Halogen bulb 12V/90W

**Light color temperature:** Up to 4,000 Kelvin

**Color:** White – RAL 9002

<b>Configurations:</b>	Ceiling	Dual C.	Wall	Mobile	Combo	FO	Battery
<b>Weight:</b>	28.0 kg	43.7kg	11.7 kg	39.2 kg	42.2 kg	0.5 kg	25.0 kg

<b>Ceiling</b>	Art. No.: <b>28506</b>
<b>Ceiling with Battery Back-up</b>	Art. No.: <b>28522</b>
<b>Dual Ceiling</b>	Art. No.: <b>28507</b>
<b>Wall</b>	Art. No.: <b>28504</b>
<b>Wall with Battery Back-up</b>	Art. No.: <b>28520</b>
<b>Mobile</b>	Art. No.: <b>28505</b>
<b>Mobile with Battery Back-up</b>	Art. No.: <b>28521</b>
<b>Combo Ceiling SunSurgery/SunSpot</b>	Art. No.: <b>28516</b>
<b>Combo Ceiling SunSurgery/SunSpot with Fiber Optic</b>	Art. No.: <b>28517</b>
<b>Single Trolley Ceiling Track</b>	Art. No.: <b>28518</b>
<b>Dual Trolley Ceiling Track</b>	Art. No.: <b>28519</b>



Dual Ceiling



Wall